

THE AIMSTRIP® PLUS BLOOD GLUCOSE METER

It provides fast, reliable blood glucose readings in 5 seconds with advanced biosensor technology, delivering precise results for confident diabetes management, and its compact, user-friendly design suits both home and professional use.



Key Features & Benefits



FAST AND ACCURATE RESULTS

Provides blood glucose readings in just 5 seconds with high accuracy.



SMALL SAMPLE SIZE

Requires only 0.7 µL of blood for testing, reducing discomfort.



NO CODING REQUIRED

Eliminates the need for manual coding, ensuring ease of use and reducing errors.



WIDE HEMATOCRIT RANGE

Works effectively with a hematocrit range of 30% to 55%.



LARGE MEMORY CAPACITY

Stores up to 250 test results with date and time for better tracking and analysis.



ALTERNATE SITE TESTING

Allows testing from fingertips, palm, or forearm for convenience and reduced pain.



AUTOMATIC SHUT-OFF

Saves battery life by turning off automatically after 2 minutes of inactivity.

+91 88078 48106

business@genevixequipments.com

www.genevixequipments.com



THE AIMSTRIP® PLUS BLOOD GLUCOSE METER



Technical Specifications

Testing Time	5 seconds
Sample Size	0.7 µL
Hematocrit Range	30% – 55%
Memory Capacity	250 results with date & time
Power Source	CR2032 3V lithium battery
Meter Size	90mm x 55mm x 15mm
Weight	47g (including battery)

Why Choose the AimStrip® Plus Blood Glucose Meter?

- ✓ Clinically tested for accuracy
- ✓ Easy to use and maintain
- ✓ Reliable performance for home and clinical settings

+91 88078 48106

business@genevixequipments.com

www.genevixequipments.com

