

ASP402 Nasal Aspirator

The ASP402 Nasal Aspirator is designed to provide quick and comfortable relief from nasal congestion for infants and young children. With its adjustable suction levels and quiet operation, it ensures a gentle yet effective nasal cleaning experience. The ergonomic design and easy-to-clean parts make it a must-have for parents looking to keep their child's nasal passages clear and healthy.



Key Features and Benefits

- ✓ **Adjustable Suction Levels**
Three suction modes to suit different levels of congestion, providing customized relief.
- ✓ **Quiet and Comfortable Operation**
Low noise design ensures a calm and stress-free experience for babies.
- ✓ **Soft Silicone Tips**
Medical-grade silicone tips are soft, gentle, and safe for delicate nasal passages.
- ✓ **Easy to Clean**
All parts are detachable and washable for hygienic maintenance.
- ✓ **USB Rechargeable**
Built-in 1200mAh battery provides long-lasting use and convenient charging.
- ✓ **One-Hand Operation**
Compact and lightweight design allows for easy handling and single-hand use.
- ✓ **LED Display**
Easy-to-read display shows the current suction level and battery status.

Technical Specifications

Suction Power	Up to 65 kPa
Suction Levels	Three adjustable levels
Battery	1200mAh rechargeable lithium battery
Noise Level	≤ 60 dB (quiet operation)
Material	Medical-grade silicone (BPA-free)
Display	LED screen
Charging	USB Type-C
Weight	135g (lightweight)



Why Choose the ASP402 Nasal Aspirator?



Clinically tested for safety and efficiency



Gentle and non-invasive



User-friendly design with easy maintenance

GET IN TOUCH



+91 88078 48106

business@genevixequipments.com

www.genevixequipments.com