

JPD-HA121

Blood Pressure Monitor

The JPD-HA121 Blood Pressure Monitor is a high-precision, fully automatic device designed for accurate and hassle-free blood pressure monitoring at home or in clinical settings. With its advanced oscillometric technology, it ensures fast and reliable measurements, making it an ideal choice for health-conscious individuals, hospitals, and medical professionals.



Key Features & Benefits



Clinically Accurate Measurements

Uses advanced oscillometric technology for precise blood pressure and pulse readings.



Large LCD Display

Clear, easy-to-read screen with bold numbers, ideal for all age groups.



One-Touch Operation

Fully automatic, enabling quick and effortless blood pressure monitoring.



Dual-User Memory Storage

Stores up to 99 readings per user, helping track health trends over time.



Irregular Heartbeat Detection

Alerts users of potential heart irregularities for early diagnosis.



Comfortable Cuff Design

Adjustable, soft cuff ensures a comfortable and Auto Shut-Off.



Why Choose the JPD-HA121 Blood Pressure Monitor?

Clinically validated for high accuracy

User-friendly, one-touch operation

Reliable readings with irregular heartbeat detection

CONTACT



+91 88078 48106



business@genevixequipments.com



www.genevixequipments.com

Technical Specifications

Measurement Method	Oscillometric method
Measurement Range	Blood Pressure: 0–280 mmHg Pulse Rate: 40–199 beats/min
Accuracy	±3 mmHg (BP), ±5% (Pulse Rate)
Display	Large LCD screen with digital readout
Memory Storage	99 readings per user (2 users)
Power Supply	4 x AA batteries or USB power
Cuff Size	22–42 cm adjustable cuff
Automatic Power Off	Yes