

FREESTYLE LIBRE 2



Key Features

Real-Time Glucose Readings

Get your current glucose reading with quick scan using your smartphone or reader. Provides glucose trends and patterns to make informed decisions.

Optional Real-Time Alerts

Customizable alarms for high or low glucose levels. Stay ahead of glucose fluctuations with proactive notifications.

Painless Sensor Application

Small, discreet sensor applied to the back of the upper arm. Comfortable to wear for up to 14 days.

No Routine Fingersticks

Accurate readings without the need for painful finger pricks. Convenient and easy to use.

Water-Resistant Sensor

Wear it while showering, swimming, or exercising. Designed for active lifestyles.

Seamless App Integration

Monitor glucose levels directly from the FreeStyle Libre app. Access glucose history, patterns, and trends anytime, anywhere.

Easy Data Sharing

Share real-time glucose data with caregivers and healthcare providers. Facilitates better diabetes management and support.



Technical Specifications

Sensor Duration	Up to 14 Days
Sensor Size	35mm x 5mm
Weight	5g
Operating Temp.	10°C – 45°C (50°F – 113°F)
Water Resistance	Up to 1 meter for 30 minutes
Measurement Range	40 to 500 mg/dL
Weight	135g (lightweight)
Connectivity	Bluetooth
Compatible Devices	iOS & Android

Why Choose FreeStyle Libre 2?

- ✓ **Effortless Monitoring**
Get real-time glucose data with a quick scan.
- ✓ **Customizable Alerts**
Be notified of high and low glucose levels instantly.
- ✓ **No Fingersticks**
Accurate readings without the discomfort of finger pricks.
- ✓ **Lifestyle Friendly**
Water-resistant sensor for everyday use.
- ✓ **Smart Insights**
Track trends and patterns directly from the app.
- ✓ **Better Health Management**
Share data with healthcare providers for enhanced care.

How It Works

Apply the Sensor

Attach the small sensor to the back of your upper arm.

Scan to Check

Use the app or reader to scan the sensor and get real-time glucose data.

Track & Manage

Monitor trends, receive alerts, and make informed health decisions.

+91 88078 48106

business@genevixequipments.com

www.genevixequipments.com

