

AimStrip® Urine Analyzer 2



Key Features & Benefits



Fast & Accurate Testing

Provides results in seconds with high accuracy



Automatic Strip Detection

Smart sensor technology ensures precise measurements



Large & Easy-to-Read LCD Display

Clear backlit screen for better visibility in any setting



User-Friendly Operation

Simple one-touch functionality for quick testing



Efficient Data Management

Stores up to 2,000 test results with date and time stamps



Compact & Portable Design

Lightweight and space-saving, perfect for small labs and clinics



Technical Specifications

Test Parameters	Compatible with multiple urine reagent strips
Throughput	Up to 120 tests per hour
Display	High-resolution LCD screen for clear result visibility
Data Storage	Stores up to 2,000 test results
Detection Method	Advanced reflectance photometry for accuracy
Connectivity	USB, LIS system integration
Power Supply	AC adapter & battery-powered options for flexibility

Why Choose AimStrip® Urine Analyzer 2?

- ✓ Accurate & high-speed testing
- ✓ User-friendly and compact design
- ✓ Seamless connectivity for better workflow

GET IN TOUCH



+91 88078 48106



business@genevixequipments.com



www.genevixequipments.com